

Graphite Crucibles and Accessories

ELTRA#	Alpha#	Description	UOM
90190	AR6247	Graphite Crucible (0.955" x 0.50")	1000
90180	AR5431P	Inner Graphite Crucible (0.623" x 0.50")	100
	AR5431P	Inner Graphite Crucible (0.623" x 0.50")	1000
90185	AR433	Outer Graphite Crucible (0.956" x 0.585")	50
88400-0491	AR053HD	Graphite Crucible (0.955" x 0.50")	1000
31360	AR31360	Graphite Electrode Tip For ELTRA OHN Analyzers	EA

Capsules, Baskets and Fluxes

ELTRA#	Alpha#	Description	UOM
90257	AR567	Smooth Wall Nickel Capsule (7 x 3.5mm) Flat Bottom	100
90256	AR569-250	Smooth Wall Nickel Capsule (10 x 5mm) Flat Bottom	250
88400-0066	AND3500	Pressed Nickel Capsule (12.5 x 5mm)	100
	AR568	Smooth Wall Nickel Capsule (14.5 x 7mm) Flat Bottom	100
	AR2183	Smooth Wall Nickel Capsule (10 x 6.5mm) Flat Bottom	100
90252	AR059	Smooth Wall Tin Capsule (H=18mm, D=5mm)	100
90250	AR2344	Nickel Basket, (1g baskets)	100
88600-0012	AR2344	Nickel Basket, (1g baskets)	100
90251	AR639	Tin Flux (0.5g pellets)	454g
90800	AR073	Graphite Powder	50g
90258	AR598	Short Nibbled Nickel Accelerator, Low Oxygen	100g


Reagents and Chemicals

ELTRA#	Alpha#	Description	UOM
90200	AR171	Alphadri® Magnesium Perchlorate - Nominal Mesh Size 10-20 (HMA HMG)	454g
90210	AR2176	AlphaSolve® II (8-20 Mesh) HAZMAT	500g
90270	AR611	Schutze Reagent (HMA HMG)	100g
90289	AR170	Rare Earth Copper Oxide	100g
90330	AR177	Quartz Wool Fine, (9 Micron)	50g
90331	AR081	Glass Wool	454g
92610	AR241	Silicone Grease	150g



Graphite Crucible (0.955" x 0.50")
AR053HD (88400-0491)

Excellent for all applications plus high melting point applications like titanium and zirconium



Smooth Wall Nickel Capsule (14.5 x 7mm) Flat Bottom
AR568



Nickel Basket, (1g baskets)
AR2344 (90250)



Tin Flux (0.5g pellets)
AR639 (90251)



Graphite Powder
AR073 (90800)



AlphaSolve® II (8-20 Mesh) HAZMAT
AR2176 (90210)



Rare Earth Copper Oxide
AR170 (90289)

Operating Supplies

ELTRA#	Alpha#	Description	UOM
88400-0472	AR138	Tweezer, Straight, Self-Closing	EA
88400-0213	AR766	Deluxe Crucible Tweezers, For Use With Graphite Crucibles	EA
71010	AR71010	ELTRA Cleaning Brush For Pedestal	EA
71032	AR938	Composite Brush	EA
71035	AR71035	ELTRA Furnace Cleaning Brush, ONH Electrode Style Furnace	EA
71031	AR14072B	ELTRA Cleaning Brush	EA

Reagent and Catalyst Tubes

ELTRA#	Alpha#	Description	UOM
88400-0006	AR11064	ELTRA Reagent Tube (280 X 20mm)	EA
11062	AR11062	ELTRA Reagent Tube (160 X 16mm)	EA
88400-0422	AR88400-0422	ELTRA Reagent Tube (240 X 20mm)	EA
20040	AR20040	ELTRA Catalyst Tube (200 X 16mm)	EA

O-rings

ELTRA#	Alpha#	Description	UOM
70230	AR70230	O-ring	EA
70210	AR70210	O-ring	EA
70435	AR70435	O-ring	EA
70430	AR70430	O-ring	EA
70405	AR70405	O-ring	EA
70410	AR70410	O-ring	EA

Oxygen, Hydrogen, and Nitrogen Standards

Material	Type	Weight	Oxygen Range	Nitrogen Range	Hydrogen Range
Steel	Pin	0.5-1.0g	3-300ppm	20-4000ppm	2-8ppm
Titanium	Pin	0.1-0.25g	800-4000ppm	50-400ppm	10-300ppm
Steel	Chip			30-1000ppm	
Zirconium	Pin	0.1g	1000-2000ppm	20-50ppm	10-300ppm
Copper	Pin	1g	8-500ppm		
Iron	Powder			0.01-3%	

please note: standard values may vary, contact Alpha or visit our website for current values

A complete list of Alpha's wide range of standards can be found on our website as well as additional helpful information, including stock status, Certificate of Analysis, and values.



Tweezer, Straight, Self-Closing
AR138 (88400-0472)



ELTRA Cleaning Brush For Pedestal
AR71010 (71010)



Composite Brush
AR938 (71032)



ELTRA Cleaning Brush
AR14072B (71031)



ELTRA Reagent Tube (160 X 16mm)
AR11062 (11062)